



Elevating Healthcare with Precision in Surgery,  
Dentistry, and Beauty Instruments for decades.

**We stand as a premier manufacturer and exporter of exceptional surgical, dental, and beauty instruments, boasting a rich legacy spanning more than two decades. Our commitment to excellence and unwavering reliability has earned us acclaim among customers and partners worldwide.**

We have a team of skilled and experienced craftsmen, engineers, and quality controllers who work together to create and deliver the best instruments for our customers. We also have a team of friendly and professional customer service representatives who are ready to assist and support our customers with their inquiries and orders.



We are committed to providing our customers with the best products and services at the most competitive prices. We also offer customized solutions and designs to meet the specific needs and preferences of our customers. We are always looking for new and innovative ways to improve our products & services and to satisfy our customers.

#### **Expertise that Matters:**


With 20 years of specialized experience, we bring unparalleled expertise to the realm of surgical instrument manufacturing. Our seasoned professionals are dedicated to delivering instruments of the highest quality.

#### **Innovation at the Core:**

Embracing cutting-edge technologies, we ensure that our manufacturing processes are at the forefront of industry innovation. This commitment guarantees that you receive instruments that meet and exceed the latest standards in medical science.


#### **Rigorous Quality Assurance:**

Quality is not just a goal; it's our standard. Our stringent quality control measures are integrated at every stage, ensuring that each instrument leaving our facility is of impeccable quality, offering reliability in every use.





## PLASTER INSTRUMENTS




**17-231**  
Nail splitting  
forceps  
135 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-231



**17-226**  
Nail extracting  
forceps  
135 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-226

aurifyinternational.com



17-225  
Cutter for cutting off  
KIRSCHNER wires,  
Tibia head screws,  
Screws, pins, up to 8 mm diam.  
470 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-225



17-223  
220 mm  
for soft wire  
up to diam. 3.5 mm  
For hard wire  
up to diam. 2.5 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-223



## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-219




## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-216

aurifyinternational.com




**17-213  
HARVEY**  
125 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-213



**17-211  
120 mm**      **17-212  
DUROTIP  
120 mm**

for soft wire  
up to diam.  
for hard wire  
up to diam.

1,2 mm  
0,8 mm

1,2 mm  
0,8 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-211

aurifyinternational.com



**17-209**  
ANGLE  
130 mm**17-210**  
DUROTIP  
ANGLE  
130 mmfor soft wire  
up to diam.  
for hard wire  
up to diam.0,5 mm  
0,3 mm0,6 mm  
0,4 mm

## PLASTER INSTRUMENTS

[View Online](#)**SKU:** AI-17-209**17-208**  
Bending pliers  
for wires up to  
1.8 mm diam.

## PLASTER INSTRUMENTS

[View Online](#)**SKU:** AI-17-208


aurifyinternational.com



## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-204



## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-203





17-200 - 17-201  
Extraction pliers  
for boring wires

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-200 - 17-201



17-197  
with serrated jaws  
17-198  
with longitudinal and  
transversal grooves  
170 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-197

aurifyinternational.com



**17-195**  
with serrated jaws  
140 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-195



**17-193**  
DUROGRIP  
Wire, tightening pliers,  
Ligature pliers  
150 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-193

aurifyinternational.com



17-239  
290 mm

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-239



17-186 - 17-188  
DEMEL  
Wire guides

## PLASTER INSTRUMENTS

[View Online](#)

**SKU:** AI-17-186 - 17-188


aurifyinternational.com



17-235  
STILLE-AESCLAP  
260 mm

PLASTER INSTRUMENTS

[View Online](#)  
SKU: AI-17-235



17-171 - 17-177  
KIRSCHNER  
Wire traction bows,  
to be used with  
KIRSCHNER-wires  
up to 2.5 mm diam.,  
With 3 hooks,  
stainless steel

|        | Inside dimensions<br>width | depth  |                       |
|--------|----------------------------|--------|-----------------------|
| 17-171 | 200 mm                     | 155 mm | for the hip           |
| 17-172 | 155 mm                     | 155 mm | for the femur         |
| 17-173 | 135 mm                     | 130 mm | for the knee          |
| 17-174 | 120 mm                     | 120 mm | for the knee          |
| 17-175 | 105 mm                     | 100 mm | for the knee          |
| 17-176 | 95 mm                      | 110 mm | for the calcaneus     |
| 17-177 | 95 mm                      | 70 mm  | for the arm and elbow |

PLASTER INSTRUMENTS

[View Online](#)  
SKU: AI-17-185

aurifyinternational.com



### Global Compliance and Recognition:

Our instruments adhere to international regulatory standards, providing you with the assurance that our products meet the highest industry benchmarks. We take pride in our global recognition as a trusted supplier in the field.



### Exceptional Value, Competitive Pricing:

We believe in offering exceptional value without compromise. Our competitive pricing ensures that you receive top-tier surgical instruments at a cost-effective rate, allowing you to optimize your budget without sacrificing quality.



### Customer-Centric Support:

Your satisfaction is our priority. Our customer support team is ready to assist you at every step – from inquiries and orders to post-purchase support. Experience a level of service that reflects our commitment to your success.



### Sustainability in Action:

Committed to environmental responsibility, Aurify embraces sustainable practices. From sourcing materials to reducing waste, we are dedicated to contributing to a greener, more sustainable future.

## We Manufacture High-Quality Instruments

### Surgical Instruments

Explore our diverse array of surgical instruments, meticulously crafted to meet a comprehensive spectrum of your medical needs.

### Beauty Instruments

Explore our extensive collection of beauty instruments crafted to meet all your aesthetic requirements.


### Dental Instruments

Broad range of dental instruments, thoughtfully designed to address a wide spectrum of oral healthcare needs. From precision drills to specialized extraction tools, our range ensures excellence in dental procedures.

### Customized Instruments

Crafting excellence to your exact specifications, our customized instruments redefine precision in every procedure. With specialized designs and personalized materials, Aurify ensures instruments that seamlessly align with your unique requirements.

### Reach Us:



15 Moosa Street, Karimpura Road  
Sialkot – 51310, Pakistan



(+92) 324 6487901



(+971) 55 531 4375



info@aurifyinternational.com

## How To Order

Choose your desired products from our website

[www.aurifyinternational.com](http://www.aurifyinternational.com)

And send inquiry using the inquiry form or the Request a Quote Form.

Email [info@aurifyinternational.com](mailto:info@aurifyinternational.com)

or call (+92) 324 6487901 with the product codes and we will provide you with a quote.

Upon approval of the quote we will process your order.